

# Find the right size to get the perfect fit with your DEHNcare personal protective equipment.

## Glove - Determine size (cm):



Please use a tape measure to measure the circumference of your hand at the knuckles (without the thumb). Your hand should be open with the fingers together.

## DEHNcare APG - normal gauntlet (cm)

## DEHNcare APG - long gauntlet

Ordering size	8	9	10	11	12	8	9	10	11	12
	M	L	XL	2XL	3XL	M	L	XL	2XL	3XL
Total length	310	315	325	330	335	440	450	455	460	465
Gauntlet length	115	115	115	115	115	250	250	250	250	250
Hand circumference	220	240	250	265	270	220	240	250	265	270